

EURON DASTO OLS

EURON
DA'S
RANGE
IS NOW
EVEN
MORE
COMPLETE.

ONYX

A LIGHTWEIGHT STREAMLINE DESIGN
WITH AN ALUMINIUM FRAME



Operator's stool with an aluminium frame and five spoke base. A 360 degree adjuster lever housed under the seat makes it easy to adapt the height in seconds. The seat is round. The anatomical curved backrest can be adjusted in depth and to three different positions (right, centre and left). Sky comes in a wide range of 24 colours, in both the thermoformed and stitched versions.

SEAT DIAMETER: 39x9 CM
BASE DIAMETER: 58 CM
HEIGHT: ADJUSTABLE 50-65 CM
WEIGHT: 11 KG

CORAL

PRACTICAL AND FUNCTIONAL WITH
AN EASY TO CLEAN POLYAMIDE FRAME



Operator's stools with a recyclable, fireproof, scratchproof and impact-proof extruded polyamide frame fitted with five self-locking wheels. The seat is made from non-deforming polyurethane that reduces the load on the back and hips. Adjustable, rounded anatomical backrest. Wide range of colours for thermoformed and stitched versions.

SEAT DIAMETER: 39x9 CM
BASE DIAMETER: 60-46 CM
HEIGHT: ADJUSTABLE 50-68 CM
WEIGHT: 13 KG

CDS301

STYLISH AND ELEGANT CHROME
PLATED STEEL FRAME



Operator's stool with a chromium-plated steel frame and five self-pivoting wheels. 360 degree height adjuster mechanism housed under the seat can be reached and operated by either hand. The anatomical curved backrest can be adjusted in depth and to three different positions. Wide range of colours for all thermoformed and stitched versions.

SEAT DIAMETER: 40x10 CM
BASE DIAMETER: 46 CM
HEIGHT: ADJUSTABLE 51-69 CM
WEIGHT: 15 KG

EURONDA

36030 Montecchio Precalcino (Vicenza) ITALY
Via dell'Artigianato, 7
tel. +39 0445 329811
fax +39 0445 865246
info@euronda.com
www.euronda.com